
() Physiomics


## Introduction


Figure 3 illustrates how the App fits within the current standard of care pathway. The therapeutic window of chemotherapy drugs is commonly established at a population level and patient dose selection is often simply scaled with Body In the current standard of care, the first docetaxel dose is selected based on patient BSA. Surface Area (BSA). Due to large inter-individual physiological variability, this A blood sample is routinely collected on the day of first injection and sometimes on the leads to a significant number of patients being under or over-dosed [1,2]. While first day of each subsequent cycle. If clinical toxicity is observed, the next dose is reduced a limited number of precision dosing techniques exist to tailor patient-specific or the treatment terminated. Docetaxel is administered in 3 week cycles until switch to treatment, they typically require costly additional tests which severely restrict another line of treatment, unacceptable toxicity or death. their use in clinical practice.

Focusing on docetaxel for advanced prostate cancer, we have developed a demonstrator for precision dosing which requires a single weekly classical blood test in the first chemotherapy cycle. It will fit within the current clinical practice to improve patient outcome at low cost.

[1] Gurney, H. (2002). British journal of cancer, 86(8), 1297 [2] Engels, F.K. et al. (2011). Clinical Cancer Research, 17(2), 353-362.







[5] McLeod, H.L., et al. (1998.). Cancer chemotherapy and pharmacology, 42(2),155-159.

## Physiomics plc<sup>©</sup> 2019

## A Precision Dosing Application For Advanced Prostate Cancer Chemotherapy

Claire Villette<sup>1</sup>, Hitesh Mistry<sup>1</sup>, Fernando Ortega<sup>1</sup>, David Orell<sup>1</sup>, Frances Brightman<sup>1</sup>, Jim Millen<sup>1</sup>, Christophe Chassagnole<sup>1</sup>.

<sup>1</sup>Physiomics plc, The Oxford Science Park, Oxford, OX4 4GA, United Kingdom. Tel. +44 (0)1865 784980. Email: <u>cchassagnole@physiomics-plc.com</u>

## Role of the App within the standard of care pathway

- Our precision dosing App requires two additional blood tests during the first chemotherapy cycle around day 7 and day 14. Just before selecting the second chemotherapy dose, patient characteristics (age, height, weight) and blood tests results are entered into the App which outputs an evaluation of neutropenia and patient median survival probability relative to the general population under the current dose (Figures 4 and 5). Simulations of dose changes can be run in seconds.
- Toxicity information provided by the App may help clinicians to identify effective, under or overdosing and modulate the dose accordingly.

- progression-free survival models, respectively.
- Antigen (PSA) as well as other biochemical markers.
- higher estimated docetaxel exposure.

					llι	IST	ra
NeutroFit	Patient	Analysis	Refine	Notes			
In the next o	cycle, the pa	atient suffered atient would s				ay(S).	
Hover over	piot to see (	coordinates					
Save a plot				Plot 2			
L Plot 3				➡ Plot 4			

# compared to population models (grey) for administered 75mg/m<sup>2</sup> dose.

х, у	In the updated cycle, the patient would suffer grade 3 neutropenia for 7.8 day Hover over plot to see coordinates x, y Save a plot		vcle, the patient would suffer grade 3 neutropenia for 7.8 day(s).
n the updated cycle, the patient would suffer grade 3 neutropenia for 7.8 day(s). Hover over plot to see coordinates x, y Bave a plot Plot 5 Plot 6	In the updated cycle, the patient would suffer grade 3 neutropenia for 7.8 day Hover over plot to see coordinates x, y Save a plot		vcle, the patient would suffer grade 3 neutropenia for 7.8 day(s).
Hover over plot to see coordinates x, y Save a plot & Plot 5 & Plot 6	Hover over plot to see coordinates x, y Save a plot		o see coordinates
x, y Save a plot Plot 5	x, y Save a plot	Plot 6	
x, y Save a plot Plot 5	x, y Save a plot	Plot 6	
Save a plot	Save a plot	Plot 6	Ł Plot 6
♣ Plot 5		Plot 6	🛃 Plot 6
♣ Plot 5		Plot 6	Le Plot 6
♣ Plot 5		Plot 6	Ł Plot 6
	₽ Plot 5	Plot 6	L Plot 6
Ł Plot 7			
ž Plot 7			
	▲ Plot 7		

Fig 5: Screenshot of the App 'dose refinement' tab. Simulation of survival probability and neutropenia in upcoming cycle for current (blue) and increased dose of 115 mg/m<sup>2</sup> (yellow). Estimated increase in median overall survival time from 440 days to 650 days.

We have developed a precision dosing App for Docetaxel in advanced prostate cancer which requires a single weekly standard blood test in the first chemotherapy cycle. Next stages include a validation clinical trial. This App has the potential, if approved by regulators, to significantly improve patient outcome and toxicity risk at low cost without disrupting the clinical treatment pathway.

## Results

Concordance levels of 0.69 and 0.6 were obtained for the overall and

Significant predictive variables include docetaxel exposure, Prostate Specific

Patients with low estimated hematologic toxicity (neutrophil count not dropping under 0.5 billion/L) presented a median overall survival time of 480 days, against 625 days for patients with higher hematologic toxicity.

Patients with low estimated docetaxel exposure (<3 µg.h/mL) presented a median overall survival time of 450 days, against 580 days for patients with

## ative underdosing case





Fig 4: Screenshot of the App 'analysis' tab. Overall and progression free survival models (top), neutropenia model for last (left) and upcoming (right) cycles. Patient models (blue) are



## Conclusion

## www.physiomics-plc.com